


ANALISIS GRANULOMETRICOS

MUESTRA: D758

UTM X: 653157


UTM Y: 3223941



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	3.87 (100 %C)	3.9%	96.1%
n°18	1 mm	13.09 (100 %C)	17.0%	83.0%
n°25	0.71 mm	9.86 (100 %C)	26.8%	73.2%
n°35	0.5 mm	8.39 (0 %C)	35.2%	64.8%
n°60	0.25 mm	23.41 (0 %C)	58.6%	41.4%
n°80	0.18 mm	12.81 (0 %C)	71.4%	28.6%
n°120	0.125 mm	18.34 (0 %C)	89.8%	10.2%
n°230	0.063 mm	8.94 (0 %C)	98.7%	1.3%
Finos	<0,063 mm	1.29 (0 %C)	100.0%	0.0%

% FINOS	1.76%
% CASCAJO	26.82%
D16	0.465
D50	0.224
D84	0.129
Media	0.273
σ	0.168
M. distribuc.	0.297

MODA ARENAS FINAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm