


ANALISIS GRANULOMETRICOS

MUESTRA: IP96

UTM X: 635646


UTM Y: 3201925



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.00 (0 %C)	0.0%	100.0%
n°18	1 mm	0.00 (0 %C)	0.0%	100.0%
n°25	0.71 mm	0.00 (0 %C)	0.0%	100.0%
n°35	0.5 mm	0.00 (0 %C)	0.0%	100.0%
n°60	0.25 mm	1.30 (0 %C)	1.3%	98.7%
n°80	0.18 mm	25.11 (0 %C)	26.4%	73.6%
n°120	0.125 mm	66.48 (0 %C)	92.9%	7.1%
n°230	0.063 mm	7.11 (0 %C)	100.0%	0.0%
Finos	<0,063 mm	0.00 (0 %C)	100.0%	0.0%

% FINOS	0.00%
% CASCAJO	0.00%
D16	0.209
D50	0.160
D84	0.132
Media	0.167
σ	0.038
M. distribuc.	0.171

MODA ARENAS FINAS



	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm