


ANALISIS GRANULOMETRICOS

MUESTRA: IP113 (REPET)

UTM X: 633006


UTM Y: 3200527



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.00 (0 %C)	0.0%	100.0%
n°18	1 mm	0.00 (0 %C)	0.0%	100.0%
n°25	0.71 mm	0.04 (50 %C)	0.0%	100.0%
n°35	0.5 mm	0.20 (30 %C)	0.2%	99.8%
n°60	0.25 mm	22.97 (0 %C)	23.2%	76.8%
n°80	0.18 mm	39.59 (0 %C)	62.8%	37.2%
n°120	0.125 mm	34.36 (0 %C)	97.2%	2.8%
n°230	0.063 mm	2.85 (0 %C)	100.0%	0.0%
Finos	<0,063 mm	0.00 (0 %C)	100.0%	0.0%

% FINOS	0.00%
% CASCAJO	0.08%
D16	0.328
D50	0.203
D84	0.146
Media	0.225
σ	0.091
M. distribuc.	0.237

MODA ARENAS FINAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm

