


## ANALISIS GRANULOMETRICOS

MUESTRA: D756

UTM X: 653454


UTM Y: 3224474



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 ( 0 %C)	0.0%	100.0%
n°10	2 mm	1.05 ( 100 %C)	1.0%	99.0%
n°18	1 mm	6.99 ( 100 %C)	8.0%	92.0%
n°25	0.71 mm	5.83 ( 100 %C)	13.9%	86.1%
n°35	0.5 mm	7.11 ( 0 %C)	21.0%	79.0%
n°60	0.25 mm	29.95 ( 0 %C)	50.9%	49.1%
n°80	0.18 mm	25.29 ( 0 %C)	76.2%	23.8%
n°120	0.125 mm	20.63 ( 0 %C)	96.9%	3.1%
n°230	0.063 mm	2.56 ( 0 %C)	99.4%	0.6%
Finos	<0,063 mm	0.58 ( 0 %C)	100.0%	0.0%

% FINOS	0.68%
% CASCAJO	13.87%
<b>D16</b>	0.444
<b>D50</b>	0.233
<b>D84</b>	0.153
<b>Media</b>	0.277
<b>σ</b>	0.145
<b>M. distribuc.</b>	0.299

### MODA ARENAS FINAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm