


ANALISIS GRANULOMETRICOS

MUESTRA: IP61


UTM X: 647254

UTM Y: 3208803



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.00 (0 %C)	0.0%	100.0%
n°18	1 mm	0.00 (0 %C)	0.0%	100.0%
n°25	0.71 mm	0.06 (50 %C)	0.1%	99.9%
n°35	0.5 mm	0.95 (0 %C)	1.0%	99.0%
n°60	0.25 mm	41.75 (0 %C)	42.8%	57.2%
n°80	0.18 mm	33.63 (0 %C)	76.4%	23.6%
n°120	0.125 mm	22.46 (0 %C)	98.9%	1.1%
n°230	0.063 mm	1.08 (0 %C)	99.9%	0.1%
Finos	<0,063 mm	0.06 (0 %C)	100.0%	0.0%

% FINOS	0.06%
% CASCAJO	0.03%
D16	0.410
D50	0.235
D84	0.161
Media	0.269
σ	0.124
M. distribuc.	0.286



MODA ARENAS FINAS

	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm