


ANALISIS GRANULOMETRICOS

MUESTRA: D729

UTM X: 647084


UTM Y: 3250414



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	18.60 (100 %C)	18.6%	81.4%
n°10	2 mm	11.20 (100 %C)	29.8%	70.2%
n°18	1 mm	11.57 (100 %C)	41.4%	58.6%
n°25	0.71 mm	5.88 (100 %C)	47.2%	52.8%
n°35	0.5 mm	5.41 (0 %C)	52.7%	47.3%
n°60	0.25 mm	13.09 (0 %C)	65.7%	34.3%
n°80	0.18 mm	10.15 (0 %C)	75.9%	24.1%
n°120	0.125 mm	16.41 (0 %C)	92.3%	7.7%
n°230	0.063 mm	7.12 (0 %C)	99.4%	0.6%
Finos	<0,063 mm	0.57 (0 %C)	100.0%	0.0%

% FINOS	1.08%
% CASCAJO	47.25%
D16	0.442
D50	0.196
D84	0.128
Media	0.255
σ	0.157
M. distribuc.	0.285

MODA ARENAS FINAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm