


## ANALISIS GRANULOMETRICOS

MUESTRA: D739


UTM X: 645539

UTM Y: 3248771



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	42.60 ( 100 %C)	42.6%	57.4%
n°10	2 mm	13.86 ( 100 %C)	56.5%	43.5%
n°18	1 mm	15.04 ( 100 %C)	71.5%	28.5%
n°25	0.71 mm	8.14 ( 100 %C)	79.6%	20.4%
n°35	0.5 mm	5.64 ( 0 %C)	85.3%	14.7%
n°60	0.25 mm	9.01 ( 0 %C)	94.3%	5.7%
n°80	0.18 mm	2.51 ( 0 %C)	96.8%	3.2%
n°120	0.125 mm	2.51 ( 0 %C)	99.3%	0.7%
n°230	0.063 mm	0.55 ( 0 %C)	99.8%	0.2%
Finos	<0,063 mm	0.16 ( 0 %C)	100.0%	0.0%

% FINOS	0.77%
% CASCAJO	79.64%
<b>D16</b>	0.589
<b>D50</b>	0.374
<b>D84</b>	0.181
<b>Media</b>	0.381
<b>σ</b>	0.204
<b>M. distribuc.</b>	0.385



## MODA ARENAS MEDIAS

	<i>D50</i>	
<i>Cantos y gravas</i>	>2	mm
<i>Arena muy gruesa</i>	1-2	mm
<i>Arena gruesa</i>	0,5-1	mm
<i>Arena media</i>	0,25-0,5	mm
<i>Arena fina</i>	0,125-0,25	mm
<i>Arena muy fina</i>	0,062-0,125	mm
<i>Fango</i>	<0,062	mm