


ANALISIS GRANULOMETRICOS

MUESTRA: D746


UTM X: 647172

UTM Y: 3244210



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	5.88 (100 %C)	5.9%	94.1%
n°18	1 mm	30.62 (100 %C)	36.5%	63.5%
n°25	0.71 mm	22.16 (100 %C)	58.7%	41.3%
n°35	0.5 mm	15.95 (0 %C)	74.6%	25.4%
n°60	0.25 mm	20.39 (0 %C)	95.0%	5.0%
n°80	0.18 mm	3.06 (0 %C)	98.1%	1.9%
n°120	0.125 mm	1.61 (0 %C)	99.7%	0.3%
n°230	0.063 mm	0.24 (0 %C)	99.9%	0.1%
Finos	<0,063 mm	0.08 (0 %C)	100.0%	0.0%

% FINOS	0.19%
% CASCAJO	58.66%
D16	0.623
D50	0.442
D84	0.270
Media	0.445
σ	0.177
M. distribuc.	0.446



MODA ARENAS MEDIAS

	<i>D50</i>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm