


ANALISIS GRANULOMETRICOS

MUESTRA: IP114

UTM X: 632436


UTM Y: 3200456



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.00 (0 %C)	0.0%	100.0%
n°18	1 mm	0.16 (50 %C)	0.2%	99.8%
n°25	0.71 mm	0.49 (50 %C)	0.7%	99.3%
n°35	0.5 mm	2.08 (25 %C)	2.7%	97.3%
n°60	0.25 mm	32.23 (0 %C)	35.0%	65.0%
n°80	0.18 mm	31.25 (0 %C)	66.2%	33.8%
n°120	0.125 mm	32.03 (0 %C)	98.2%	1.8%
n°230	0.063 mm	1.72 (0 %C)	100.0%	0.0%
Finos	<0,063 mm	0.04 (0 %C)	100.0%	0.0%

% FINOS	0.04%
% CASCAJO	0.85%
D16	0.392
D50	0.215
D84	0.149
Media	0.252
σ	0.121
M. distribuc.	0.270

MODA
ARENAS FINAS



	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm