


## ANALISIS GRANULOMETRICOS

MUESTRA: D732

UTM X: 647184


UTM Y: 3251389



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 ( 0 %C)	0.0%	100.0%
n°10	2 mm	0.20 ( 100 %C)	0.2%	99.8%
n°18	1 mm	1.30 ( 100 %C)	1.5%	98.5%
n°25	0.71 mm	2.21 ( 100 %C)	3.7%	96.3%
n°35	0.5 mm	4.21 ( 0 %C)	7.9%	92.1%
n°60	0.25 mm	50.05 ( 0 %C)	58.0%	42.0%
n°80	0.18 mm	23.97 ( 0 %C)	81.9%	18.1%
n°120	0.125 mm	13.84 ( 0 %C)	95.8%	4.2%
n°230	0.063 mm	3.41 ( 0 %C)	99.2%	0.8%
Finos	<0,063 mm	0.80 ( 0 %C)	100.0%	0.0%

% FINOS	0.83%
% CASCAJO	3.71%
<b>D16</b>	0.444
<b>D50</b>	0.281
<b>D84</b>	0.169
<b>Media</b>	0.298
<b>σ</b>	0.137
<b>M. distribuc.</b>	0.307

## MODA ARENAS MEDIAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm