


ANALISIS GRANULOMETRICOS

MUESTRA: IP083


UTM X: 647895

UTM Y: 3209198



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.28 (100 %C)	0.3%	99.7%
n°18	1 mm	0.09 (100 %C)	0.4%	99.6%
n°25	0.71 mm	0.09 (100 %C)	0.5%	99.5%
n°35	0.5 mm	0.57 (0 %C)	1.0%	99.0%
n°60	0.25 mm	23.25 (0 %C)	24.3%	75.7%
n°80	0.18 mm	38.75 (0 %C)	63.0%	37.0%
n°120	0.125 mm	31.66 (0 %C)	94.7%	5.3%
n°230	0.063 mm	5.20 (0 %C)	99.9%	0.1%
Finos	<0,063 mm	0.09 (0 %C)	100.0%	0.0%

% FINOS	0.09%
% CASCAJO	0.47%
D16	0.335
D50	0.203
D84	0.143
Media	0.227
σ	0.096
M. distribuc.	0.239



MODA ARENAS FINAS

	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm