


ANALISIS GRANULOMETRICOS

MUESTRA: D742


UTM X: 644018

UTM Y: 3244338



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	7.75 (100 %C)	7.8%	92.2%
n°10	2 mm	8.66 (100 %C)	16.4%	83.6%
n°18	1 mm	34.24 (100 %C)	50.6%	49.4%
n°25	0.71 mm	19.90 (100 %C)	70.5%	29.5%
n°35	0.5 mm	12.14 (0 %C)	82.7%	17.3%
n°60	0.25 mm	13.05 (0 %C)	95.7%	4.3%
n°80	0.18 mm	1.94 (0 %C)	97.7%	2.3%
n°120	0.125 mm	1.42 (0 %C)	99.1%	0.9%
n°230	0.063 mm	0.52 (0 %C)	99.6%	0.4%
Finos	<0,063 mm	0.39 (0 %C)	100.0%	0.0%

% FINOS	1.32%
% CASCAJO	70.54%
D16	0.629
D50	0.450
D84	0.259
Media	0.446
σ	0.185
M. distribuc.	0.444



MODA ARENAS MEDIAS

	<i>D50</i>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm