


ANALISIS GRANULOMETRICOS

MUESTRA: IP106

UTM X: 633709


UTM Y: 3200675



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.00 (0 %C)	0.0%	100.0%
n°18	1 mm	0.04 (50 %C)	0.0%	100.0%
n°25	0.71 mm	0.09 (50 %C)	0.1%	99.9%
n°35	0.5 mm	0.26 (50 %C)	0.4%	99.6%
n°60	0.25 mm	14.81 (0 %C)	15.2%	84.8%
n°80	0.18 mm	43.38 (0 %C)	58.6%	41.4%
n°120	0.125 mm	38.62 (0 %C)	97.2%	2.8%
n°230	0.063 mm	2.80 (0 %C)	100.0%	0.0%
Finos	<0,063 mm	0.00 (0 %C)	100.0%	0.0%

% FINOS	0.00%
% CASCAJO	0.20%
D16	0.248
D50	0.194
D84	0.144
Media	0.195
σ	0.052
M. distribuc.	0.196

MODA ARENAS FINAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm