


ANALISIS GRANULOMETRICOS

MUESTRA: D736

UTM X: 646829


UTM Y: 3249949



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	14.36 (100 %C)	14.4%	85.6%
n°10	2 mm	4.75 (100 %C)	19.1%	80.9%
n°18	1 mm	5.52 (100 %C)	24.6%	75.4%
n°25	0.71 mm	3.54 (100 %C)	28.2%	71.8%
n°35	0.5 mm	3.09 (0 %C)	31.3%	68.7%
n°60	0.25 mm	20.99 (0 %C)	52.3%	47.7%
n°80	0.18 mm	20.11 (0 %C)	72.4%	27.6%
n°120	0.125 mm	20.55 (0 %C)	92.9%	7.1%
n°230	0.063 mm	5.75 (0 %C)	98.7%	1.3%
Finos	<0,063 mm	1.33 (0 %C)	100.0%	0.0%

% FINOS	1.85%
% CASCAJO	28.18%
D16	0.400
D50	0.209
D84	0.137
Media	0.249
σ	0.132
M. distribuc.	0.268

MODA ARENAS FINAS



	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm