


ANALISIS GRANULOMETRICOS

MUESTRA: D747

UTM X: 646414


UTM Y: 3243343



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.76 (100 %C)	0.8%	99.2%
n°10	2 mm	6.61 (100 %C)	7.4%	92.6%
n°18	1 mm	28.05 (100 %C)	35.4%	64.6%
n°25	0.71 mm	20.59 (100 %C)	56.0%	44.0%
n°35	0.5 mm	16.05 (0 %C)	72.0%	28.0%
n°60	0.25 mm	21.15 (0 %C)	93.2%	6.8%
n°80	0.18 mm	4.15 (0 %C)	97.4%	2.6%
n°120	0.125 mm	2.27 (0 %C)	99.6%	0.4%
n°230	0.063 mm	0.28 (0 %C)	99.9%	0.1%
Finos	<0,063 mm	0.09 (0 %C)	100.0%	0.0%

% FINOS	0.21%
% CASCAJO	56.00%
D16	0.618
D50	0.430
D84	0.253
Media	0.433
σ	0.183
M. distribuc.	0.435

MODA ARENAS MEDIAS



	<i>D50</i>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm